

Product Datasheet

IgG antibody (HRP) (orb436182)

Description	Rabbit polyclonal antibody against IgG conjugated to HRP
Species/Host	Rabbit
Reactivity	Goat
Conjugation	HRP
Tested Applications	ELISA, IH-Fr, WB
Immunogen	Goat IgG, whole molecule
Target	IgG
Preservatives	0.01%Gentamicin sulphate1%Bovine Serum Albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Horseradish Peroxidase (HRP) - lyophilized
Concentration	2.0 mg/ml
Storage	-20°C after reconstitution
Note	For research use only
Application notes	(H/L)
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.