


Product Datasheet

HLA-G antibody (PE) (orb43619)

Description	Mouse Monoclonal to HLA-G.
Reactivity	Human
Conjugation	PE
Tested Applications	FC
Immunogen	C-terminal amino acid sequence (22-mer) of soluble HLA-G5 and HLA-G6 proteins coupled to ovalbumin.
Target	HLA-G
Preservatives	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	0.1 mg/ml
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 2-3 µg/ml.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified antibody is conjugated with R-phycoerythrin (PE) under optimum conditions. Unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.
Uniprot ID	P17693
Entrez	3135
Dilution Range	Flow cytometry: Recommended dilution: 1-4 µg/ml. Intracellular staining.
Expiration Date	12 months from date of receipt.



Separation of
HLA-G
transfected
HEK-293 ...