

Product Datasheet

IgG antibody (HRP) (orb436204)

Description	Goat polyclonal antibody against IgG conjugated to HRP
Species/Host	Goat
Reactivity	Rabbit
Conjugation	HRP
Tested Applications	ELISA, IH-Fr, WB
Immunogen	Native rabbit IgG.
Target	IgG
Preservatives	0.01% Gentamicin sulphate10 mg/ml Bovine Serum Albumin. Phosphate buffered saline
Form/Appearance	Ig Fraction conjugated to Horseradish Peroxidase (HRP) - lyophilized
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	(Fc)
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.