

Product Datasheet

IgM antibody (HRP) (orb436216)

Description	Goat polyclonal antibody against IgM conjugated to HRP
Species/Host	Goat
Reactivity	Mouse
Conjugation	HRP
Tested Applications	ELISA, IH-Fr, WB
Immunogen	Mouse IgM mu chain.
Target	IgM
Preservatives	0.01%Gentamicin sulphate1%Bovine Serum Albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Horseradish Peroxidase (HRP) - lyophilized
Storage	-20°C after reconstitution
Note	For research use only
Application notes	MU CHAIN
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.