



Product Datasheet

Blood Group Lewis a antibody (orb43622)

Description	Mouse Monoclonal to Blood Group Lewis a.
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	Mucin isolated by size-exclusion chromatography from an ovarian cyst fluid (pure endocervical type according to the Fenoglio's classification).
Target	Blood Group Lewis a
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Immunohistochemistry (paraffin sections): Recommended dilution: 5 µg/ml; positive control: colorectal adenocarcinoma; heat mediated antigen retrieval (sodium citrate); mAb incubation 1 hour / RT, detection DAB.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Dilution Range	Immunohistochemistry (paraffin sections): Recommended dilution: 5 µg/ml; positive control: colorectal adenocarcinoma; heat mediated antigen retrieval (sodium citrate); mAb incubation 1 hour / RT, detection DAB.
Expiration Date	12 months from date of receipt.



Immunohistochemistry staining of human c...



Immunohistochemistry staining of human p...