

Product Datasheet

IgG1 antibody (AP) (orb436228)

Description	Sheep polyclonal antibody against IgG1 conjugated to AP
Species/Host	Sheep
Reactivity	Bovine
Conjugation	AP
Tested Applications	ELISA
Immunogen	Purified bovine IgG1.
Target	IgG1
Preservatives	0.09% Sodium Azide (NaN ₃). 50mM HEPES, 0.1M NaCl, 1mM MgCl ₂ , 0.1mM ZnCl ₂
Form/Appearance	Purified IgG conjugated to Alkaline Phosphatase - liquid
Concentration	0.5mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.