

Product Datasheet

IgG2 antibody (FITC) (orb436234)

Description	Sheep polyclonal antibody against IgG2 conjugated to FITC
Species/Host	Sheep
Reactivity	Bovine
Conjugation	FITC
Tested Applications	FC, IH-Fr
Immunogen	Purified bovine IgG2.
Target	IgG2
Preservatives	0.09% Sodium Azide. Phosphate buffered saline
Form/Appearance	Purified IgG fraction conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Concentration	1.0 mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.