

Product Datasheet

IgG antibody (orb436236)

Description	Sheep polyclonal antibody to IgG
Species/Host	Sheep
Reactivity	Bovine
Conjugation	Unconjugated
Tested Applications	ELISA, ID
Immunogen	Purified bovine IgG.
Target	IgG
Preservatives	0.09% Sodium Azide. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.