

## Product Datasheet

### IgA antibody (HRP) (orb436245)

|                            |                                                                         |
|----------------------------|-------------------------------------------------------------------------|
| <b>Description</b>         | Goat polyclonal antibody against IgA conjugated to HRP                  |
| <b>Species/Host</b>        | Goat                                                                    |
| <b>Reactivity</b>          | Feline                                                                  |
| <b>Conjugation</b>         | HRP                                                                     |
| <b>Tested Applications</b> | ELISA, WB                                                               |
| <b>Immunogen</b>           | Purified feline IgA.                                                    |
| <b>Target</b>              | IgA                                                                     |
| <b>Preservatives</b>       | 0.05% Proclin 300 0.2% Bovine Serum Albumin. Phosphate buffered saline. |
| <b>Form/Appearance</b>     | Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid        |
| <b>Concentration</b>       | 1.0 mg/ml                                                               |
| <b>Storage</b>             | +4°C DO NOT FREEZE                                                      |
| <b>Note</b>                | For research use only                                                   |
| <b>Isotype</b>             | Polyclonal IgG                                                          |
| <b>Clonality</b>           | Polyclonal                                                              |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                         |