

Product Datasheet

IgM antibody (FITC) (orb436250)

Description	Goat polyclonal antibody against IgM conjugated to FITC
Species/Host	Goat
Reactivity	Gallus
Conjugation	FITC
Tested Applications	FC, IH-Fr
Immunogen	Purified chicken IgM.
Target	IgM
Preservatives	0.09% Sodium Azide 0.2% Bovine Serum Albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Concentration	1.0 mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.