

Product Datasheet

IgA antibody (FITC) (orb436254)

Description	Goat polyclonal antibody against IgA conjugated to FITC
Species/Host	Goat
Reactivity	Gallus
Conjugation	FITC
Tested Applications	FC, IH-Fr
Immunogen	Purified chicken IgA.
Target	IgA
Preservatives	0.09% Sodium Azide 0.2% Bovine Serum Albumin. Phosphate buffered saline
Form/Appearance	Purified IgG fraction conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Concentration	1.0 mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.