

Product Datasheet

IgG antibody (orb436256)

Description	Goat polyclonal antibody to IgG
Species/Host	Goat
Reactivity	Gallus
Conjugation	Unconjugated
Tested Applications	ELISA, ID
Immunogen	Purified chicken IgG (Fc) fragment.
Target	IgG
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Application notes	(Fc)
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.