

Product Datasheet

IgG antibody (AP) (orb436257)

Description	Goat polyclonal antibody against IgG conjugated to AP
Species/Host	Goat
Reactivity	Gallus
Conjugation	AP
Tested Applications	ELISA
Immunogen	Purified chicken IgG (Fc) fragment.
Target	IgG
Preservatives	0.09% Sodium Azide 0.2% Bovine Serum Albumin. 50mM HEPES, 0.1M NaCl, 1mM MgCl ₂ , 0.1mM ZnCl ₂
Form/Appearance	Purified IgG conjugated to Alkaline Phosphatase - liquid
Concentration	0.5mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Application notes	(Fc)
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.