

Product Datasheet

IgG antibody (FITC) (orb436258)

Description	Goat polyclonal antibody against IgG conjugated to FITC
Species/Host	Goat
Reactivity	Gallus
Conjugation	FITC
Tested Applications	FC
Immunogen	Purified chicken IgG (Fc) fragment.
Target	IgG
Preservatives	0.09% Sodium Azide (NaN ₃) 0.2% Bovine Serum Albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Concentration	1.0 mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Application notes	(Fc)
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.