

Product Datasheet

IgM antibody (orb436260)

Description	Goat polyclonal antibody to IgM
Species/Host	Goat
Reactivity	Canine
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IH-Fr
Immunogen	Purified canine IgM.
Target	IgM
Preservatives	0.09% Sodium Azide. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.