

Product Datasheet

IgM antibody (AP) (orb436261)

Description	Goat polyclonal antibody against IgM conjugated to AP
Species/Host	Goat
Reactivity	Canine
Conjugation	AP
Tested Applications	ELISA
Immunogen	Purified canine IgM.
Target	IgM
Preservatives	0.09% Sodium Azide 0.2% Bovine Serum Albumin. 50mM HEPES, 0.1M NaCl, 1mM MgCl ₂ , 0.1mM ZnCl ₂
Form/Appearance	Purified IgG conjugated to Alkaline Phosphatase - liquid
Concentration	0.5mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.