

Product Datasheet

IgA antibody (orb436264)

Description	Goat polyclonal antibody to IgA
Species/Host	Goat
Reactivity	Canine
Conjugation	Unconjugated
Tested Applications	ELISA, ID, IH-Fr, WB
Immunogen	Purified canine IgA.
Target	IgA
Preservatives	0.09% Sodium Azide. Phosphate buffered saline
Form/Appearance	Purified IgG - Liquid
Concentration	1.0 mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.