

Product Datasheet

IgG antibody (Biotin) (orb436287)

Description	Goat polyclonal antibody against IgG conjugated to Biotin
Species/Host	Goat
Reactivity	Porcine
Conjugation	Biotin
Tested Applications	ELISA, WB
Immunogen	Purified porcine IgG.
Target	IgG
Preservatives	0.09% Sodium Azide (NaN ₃) 0.2% Bovine Serum Albumin. Phosphate buffered saline
Form/Appearance	Purified IgG fraction conjugated to Biotin - liquid
Concentration	1.0 mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Application notes	(Fc)
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.