

Product Datasheet

IgG antibody (HRP) (orb436289)

Description	Goat polyclonal antibody against IgG conjugated to HRP
Species/Host	Goat
Reactivity	Porcine
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Purified porcine IgG.
Target	IgG
Preservatives	0.05% Pro-Clean 400 0.2% Bovine Serum Albumin. Phosphate buffered saline.
Form/Appearance	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
Concentration	1.0 mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Application notes	(Fc)
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.