

## Product Datasheet

### IgG antibody (FITC) (orb436292)

|                            |                                                                                |
|----------------------------|--------------------------------------------------------------------------------|
| <b>Description</b>         | Goat polyclonal antibody against IgG conjugated to FITC                        |
| <b>Species/Host</b>        | Goat                                                                           |
| <b>Reactivity</b>          | Monkey                                                                         |
| <b>Conjugation</b>         | FITC                                                                           |
| <b>Tested Applications</b> | FC, IF, IH-Fr                                                                  |
| <b>Immunogen</b>           | Purified monkey IgG.                                                           |
| <b>Target</b>              | IgG                                                                            |
| <b>Preservatives</b>       | 0.09% Sodium Azide 0.2% Bovine Serum Albumin. Phosphate buffered saline        |
| <b>Form/Appearance</b>     | Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid |
| <b>Concentration</b>       | 1.0 mg/ml                                                                      |
| <b>Storage</b>             | +4°C DO NOT FREEZE                                                             |
| <b>Note</b>                | For research use only                                                          |
| <b>Application notes</b>   | (H/L)                                                                          |
| <b>Isotype</b>             | Polyclonal IgG                                                                 |
| <b>Clonality</b>           | Polyclonal                                                                     |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                |