

Product Datasheet

IgG antibody (orb436294)

Description	Goat polyclonal antibody to IgG
Species/Host	Goat
Reactivity	Equine
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Purified equine IgG.
Target	IgG
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Application notes	(H/L)
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.