

Product Datasheet

IgM antibody (HRP) (orb436301)

Description	Rabbit polyclonal antibody against IgM conjugated to HRP
Species/Host	Rabbit
Reactivity	Goat
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Purified goat IgM.
Target	IgM
Preservatives	0.2% Bovine Serum Albumin. 0.05% Pro-clean 400. Phosphate buffered saline.
Form/Appearance	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
Concentration	1.0 mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.