

Product Datasheet

IgA antibody (AP) (orb436312)

Description	Sheep polyclonal antibody against IgA conjugated to AP
Species/Host	Sheep
Reactivity	Bovine
Conjugation	AP
Tested Applications	ELISA, ICC, IF, IHC-Fr, WB
Immunogen	Purified bovine IgA
Target	IgA
Preservatives	0.09% Sodium Azide (NaN ₃) 0.2% Bovine Serum Albumin. 50mM HEPES, 0.1M NaCl, 1mM MgCl ₂ , 0.1mM ZnCl ₂
Form/Appearance	Purified IgG conjugated to Alkaline Phosphatase - liquid
Concentration	0.5 mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.