

Product Datasheet

IgA antibody (Biotin) (orb436313)

Description	Sheep polyclonal antibody against IgA conjugated to Biotin
Species/Host	Sheep
Reactivity	Bovine
Conjugation	Biotin
Tested Applications	ELISA, ICC, IF, IHC-Fr, WB
Immunogen	Purified bovine IgA
Target	IgA
Preservatives	0.09% Sodium Azide (NaN ₃) 0.2% Bovine Serum Albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Biotin - liquid
Concentration	1.0 mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.