

Product Datasheet

IgA antibody (FITC) (orb436314)

Description	Sheep polyclonal antibody against IgA conjugated to FITC
Species/Host	Sheep
Reactivity	Bovine
Conjugation	FITC
Tested Applications	ICC, IF
Immunogen	Purified bovine IgA
Target	IgA
Preservatives	0.09% Sodium Azide (NaN ₃) 0.2% Bovine Serum Albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Concentration	1.0 mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.