

Product Datasheet

HISTIDINE TAG antibody (FITC) (orb436319)

Description	Goat polyclonal antibody against HISTIDINE TAG conjugated to FITC
Species/Host	Goat
Conjugation	FITC
Tested Applications	ICC, IF
Immunogen	Keyhole limpet haemocyanin (KLH) conjugated H-H-H-H-H
Target	HISTIDINE TAG
Preservatives	0.09% Sodium Azide (NaN ₃) 0.2% Bovine Serum Albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Concentration	1.0mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.