

Product Datasheet

IgG antibody (HRP) (orb436322)

Description	Rabbit polyclonal antibody against IgG conjugated to HRP
Species/Host	Rabbit
Reactivity	Bovine
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Purified intact bovine IgG
Target	IgG
Preservatives	0.05% Proclin 300.. Phosphate buffered saline.
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Application notes	(H/L)
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.