

Product Datasheet

IgE antibody (orb436329)

Description	Goat polyclonal antibody to IgE
Species/Host	Goat
Reactivity	Canine
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Native canine IgE.
Target	IgE
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.