

Product Datasheet

IgG2 antibody (HRP) (orb436338)

Description	Sheep polyclonal antibody against IgG2 conjugated to HRP
Species/Host	Sheep
Reactivity	Canine
Conjugation	HRP
Tested Applications	ELISA, IH-Fr, WB
Target	IgG2
Preservatives	0.05% Proclin 300 0.2% Bovine Serum Albumin. Phosphate buffered saline.
Form/Appearance	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
Concentration	1.0 mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.