

Product Datasheet

IgG2b antibody (Biotin) (orb436349)

Description	Mouse monoclonal antibody against IgG2b conjugated to Biotin
Species/Host	Mouse
Reactivity	Rat
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Rat IgG2b.
Target	IgG2b
Preservatives	0.09% sodium azide (NaN ₃)50% glycerol. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to biotin - liquid
Concentration	1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	HEAVY CHAIN
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.