

Product Datasheet

KAPPA/LAMBDA LIGHT CHAIN antibody (HRP) (orb436353)

Description	Mouse monoclonal antibody against KAPPA/LAMBDA LIGHT CHAIN conjugated to HRP
Species/Host	Mouse
Reactivity	Rat
Conjugation	HRP
Tested Applications	ELISA
Immunogen	MARK-1: Rat kappa light chains (IR202 and IR968) MARL-15: Rat lambda light chains (RH58, IR31)
Target	KAPPA/LAMBDA LIGHT CHAIN
Preservatives	0.01% Thiomersal50% glycerol. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
Concentration	1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.