

Product Datasheet

IgG2a antibody (orb436358)

Description	Rat monoclonal antibody to IgG2a
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Purified IgG2a from Balb/c mice. PTR
Target	IgG2a
Preservatives	0.1% sodium azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	HEAVY CHAIN
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.