

## Product Datasheet

### IgG2a antibody (HRP) (orb436361)

|                            |                                                                  |
|----------------------------|------------------------------------------------------------------|
| <b>Description</b>         | Rat monoclonal antibody against IgG2a conjugated to HRP          |
| <b>Species/Host</b>        | Rat                                                              |
| <b>Reactivity</b>          | Mouse                                                            |
| <b>Conjugation</b>         | HRP                                                              |
| <b>Tested Applications</b> | ELISA                                                            |
| <b>Immunogen</b>           | Purified IgG2a from Balb/c mice. PTR                             |
| <b>Target</b>              | IgG2a                                                            |
| <b>Preservatives</b>       | 0.01% Thiomersal50% glycerol. Phosphate buffered saline          |
| <b>Form/Appearance</b>     | Purified IgG conjugated to horseradish peroxidase (HRP) - liquid |
| <b>Concentration</b>       | 1.0 mg/ml                                                        |
| <b>Storage</b>             | +4°C, -20°C if preferred                                         |
| <b>Note</b>                | For research use only                                            |
| <b>Application notes</b>   | HEAVY CHAIN                                                      |
| <b>Isotype</b>             | IgG1                                                             |
| <b>Clonality</b>           | Monoclonal                                                       |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                  |