

Product Datasheet

IgG2b antibody (orb436362)

Description	Rat monoclonal antibody to IgG2b
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Polyclonal IgG from BALB/c mice.
Target	IgG2b
Preservatives	0.09% sodium azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	HEAVY CHAIN
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.