

Product Datasheet

IgM antibody (orb436365)

Description	Mouse monoclonal antibody to IgM
Species/Host	Mouse
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	IR202 rat IgM myeloma protein
Target	IgM
Preservatives	0.1% Sodium Azide. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	HEAVY CHAIN
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.