

Product Datasheet

IgE antibody (orb436368)

Description	Mouse monoclonal antibody to IgE
Species/Host	Mouse
Reactivity	Canine
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Affinity purified IgE preparation from heavily parasitized and allergic dog serum.
Target	IgE
Preservatives	0.1% sodium azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.