

Product Datasheet

IgM antibody (FITC) (orb436372)

Description	Mouse monoclonal antibody against IgM conjugated to FITC
Species/Host	Mouse
Reactivity	Rat
Conjugation	FITC
Tested Applications	FC
Immunogen	IR202 rat IgM myeloma protein
Target	IgM
Preservatives	0.1% Sodium Azide 50% Glycerol. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	HEAVY CHAIN
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.