

Product Datasheet

IgGb antibody (orb436375)

Description Mouse monoclonal antibody to IgGb

Species/Host Mouse

Reactivity Equine

Conjugation Unconjugated

Tested Applications ELISA

Immunogen Purified equine IgGb.

Target IgGb

Preservatives 0.09% sodium azide (NaN₃)

Form/Appearance Liquid: Phosphate buffered saline

Concentration 1.0 mg/ml

Storage This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C. Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Note For research use only


Isotype IgG1

Clonality Monoclonal

Purity Purified IgG prepared by affinity chromatography on Protein A from tissue culture supernatant

Dilution Range ELISA: 1/5000-1/50000

Expiration Date 12 months from date of receipt.



Purified Horse IgG detected with Mouse a...