

Product Datasheet

IgE antibody (Biotin) (orb436382)

Description	Mouse monoclonal antibody against IgE conjugated to Biotin
Species/Host	Mouse
Reactivity	Rat
Conjugation	Biotin
Tested Applications	ELISA, WB
Immunogen	Rat IR162, IR1016, IR2 and IR410 IgE myeloma proteins.
Target	IgE
Preservatives	0.1% Sodium Azide 50% Glycerol. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Biotin - liquid
Concentration	1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.