

Product Datasheet

IgM antibody (orb436391)

Description	Rat monoclonal antibody to IgM
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IH-Fr
Immunogen	Polyclonal and monoclonal mouse IgM.
Target	IgM
Preservatives	0.1% Sodium Azide. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	HEAVY CHAIN
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.