

Product Datasheet

IgM antibody (FITC) (orb436393)

Description	Rat monoclonal antibody against IgM conjugated to FITC
Species/Host	Rat
Reactivity	Mouse
Conjugation	FITC
Tested Applications	FC
Immunogen	Polyclonal and monoclonal mouse IgM.
Target	IgM
Preservatives	0.1% Sodium Azide 50% Glycerol. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Concentration	1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	HEAVY CHAIN
Isotype	IgG2a
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.