

Product Datasheet

IgG4 antibody (orb436394)

Description	Mouse monoclonal antibody to IgG4
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	DOT, IHC-Fr, IHC-P, WB
Immunogen	Purified IgG4
Target	IgG4
Preservatives	0.1% sodium azide (NaN ₃). Borate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	0.5mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.