

Product Datasheet

LAMBDA LIGHT CHAIN antibody (orb436404)

Description	Mouse monoclonal antibody to LAMBDA LIGHT CHAIN
Species/Host	Mouse
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IH-Fr, WB
Immunogen	Rat lambda light chain.
Target	LAMBDA LIGHT CHAIN
Preservatives	0.09% sodium azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.