

Product Datasheet

IgA antibody (orb436406)

Description	Mouse monoclonal antibody to IgA
Species/Host	Mouse
Reactivity	Monkey
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Purified IgA from serum.
Target	IgA
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.