

Product Datasheet

IgA antibody (Biotin) (orb436407)

Description	Mouse monoclonal antibody against IgA conjugated to Biotin
Species/Host	Mouse
Reactivity	Monkey
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Purified IgA from serum.
Target	IgA
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Biotin - liquid
Concentration	0.5mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.