

Product Datasheet

IgG2a antibody (FITC) (orb436422)

Description	Rat monoclonal antibody against IgG2a conjugated to FITC
Species/Host	Rat
Reactivity	Mouse
Conjugation	FITC
Tested Applications	FC
Immunogen	Purified IgG from BALB/c mice.
Target	IgG2a
Preservatives	0.1% Sodium Azide 50% Glycerol. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	HEAVY CHAIN
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.