

Product Datasheet

IgG2a antibody (HRP) (orb436423)

Description	Rat monoclonal antibody against IgG2a conjugated to HRP
Species/Host	Rat
Reactivity	Mouse
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Purified IgG from BALB/c mice.
Target	IgG2a
Preservatives	0.01%Thiomersal50%Glycerol. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
Concentration	1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	HEAVY CHAIN
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.