

Product Datasheet

KAPPA LIGHT CHAIN antibody (HRP) (orb436438)

Description	Mouse monoclonal antibody against KAPPA LIGHT CHAIN conjugated to HRP
Species/Host	Mouse
Reactivity	Rat
Conjugation	HRP
Tested Applications	ELISA
Immunogen	F(ab') ₂ fragment of rat IgG.
Target	KAPPA LIGHT CHAIN
Preservatives	0.01%Thiomersal. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	(IgKa PREFERENCE)
Isotype	IgG2a
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.