

Product Datasheet

IgA1 antibody (orb436448)

Description	Mouse monoclonal antibody to IgA1
Species/Host	Mouse
Reactivity	Duck
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, WB
Immunogen	Purified duck IgA.
Target	IgA1
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	HEAVY CHAIN
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.