

Product Datasheet

IgG antibody (HRP) (orb436467)

Description	Mouse monoclonal antibody against IgG conjugated to HRP
Species/Host	Mouse
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA, IH-Fr
Immunogen	Human IgG Polyclonal.
Target	IgG
Preservatives	0.01% Thiomersal. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	(Fc) (CH2 DOMAIN)
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.